

Certificate of Analysis

Product Name	CELLPACK
Product Code	PK-30L
P. S. Code	P88408711
Lot No.	AM0020
Exp. Date (YYYY-MM-DD)	2021/09/20
Manufacturing Date (YYYY-MM-DD)	2020/03/20

Item of analysis and testing	Acceptance criterion	Result	Pass / Fail
pH	7.75 - 7.85	7.78	Pass
Electric Conductivity	13.20 - 13.40 (mS/cm)	13.26	Pass
Background Count value	$Bx \leq 110$	0	Pass

This is to certify that above product has been produced and passed the inspection according to the specification of SYSMEX Asia Pacific.



Jesslyn Tuak

Executive, Quality Control & HSE

2020/03/20